

HYDRION® PLASTIC pH STRIPS

Hydrion brand Plastic pH Strips provide instant, easy-to-read results, with clear, bright single color matches across the entire pH range.

Hydrion Plastic pH Strips are packaged in convenient flip-top vials, with matching color charts. Strips are available in 5 pH ranges. Short-range strips offer clear, bright, single-color matches at every 0.2/0.3 interval (5.5-8.0) or every 0.5 interval (0.0-6.0, 5.0-9.0, 6.5-13.0, 3.0-5.5, 3.5-6.5). The wide-range strip provides accurate pH measurement at each 1.0 interval over the entire pH range (0.0-14.0).

- Fast, accurate, easy-to-use: one dip-one color match
- Packaged in convenient flip-top vials
- Made in the U.S.A.



For more information contact Rob Lynch at 718.338.3618 Rob@microessentiallab.com www.MicroEssentialLab.com



HYDRION® PLASTIC pH STRIPS

- Strip dimensions: 2 3/4 in. x 1/4 in.
- Packaged in flip-top vials
- 100 strips per vial

Catalog No.	pH Range	pH Readings
9200	0.0-6.0	0.0 0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0 5.5 6.0
9400	5.0-9.0	5.0 5.5 6.0 6.5 7.0 7.5 8.0 8.5 9.0
9700	5.5-8.0	5.5 5.8 6.0 6.2 6.4 6.6 7.0 7.2 7.4 7.6 8.0
9600	6.5-13.0	6.5 7.0 7.5 8.0 8.5 9.0 9.5 10.0 10.5 11.0 11.5 12.0 12.5 13.0
9800	0.0-14.0	0.0 1.0 2.0 3.0 4.0 5.0 6.0 7.0 8.0 9.0 10.0 11.0 12.0 13.0 14.0
93565	3.5-6.5	3.5 4.0 4.5 5.0 5.5 6.0 6.5
93055	3.0-5.5	3.0 3.5 4.0 4.5 5.0 5.5
Sale Unit	Price	Minimum Order Quantity
Case of 6 vials	Distributor p	pricing available 1 case

How to Measure pH with Hydrion Strips

- Dip the strip into the test solution
- Remove the strip and immediately compare the resulting color with the matching pH color chart



© 2022 Micro Essential Laboratory Inc. All rights reserved. Printed in U.S.A.

For more information contact Rob Lynch at 718.338.3618 Rob@microessentiallab.com www.MicroEssentialLab.com ME9000/04.11.22