

Hydrion® Peracetic Acid (PAA) Test Strips, 0-160 and 0-2000 PPM

A Quick, reliable, and easy way to test the strength and efficacy of PAA solutions used in many applications including healthcare and food service.

PAA-2000 measures the concentration of PAA at the ranges commonly used to disinfect in hospitals, nursing homes, and other healthcare facilities.

PAA-160 measures the concentration of PAA at the common ranges used for sanitization in vegetable and fruit wash applications in food service establishments.

- A quick, simple, reliable way to ensure compliance with industry specific guidelines.
- Hydrion products are trusted by health agencies around the world
- Made in USA



PAA-160



PAA-2000



PAA-1000

For more information contact
Rob Lynch at 718.338.3618
Rob@microessentiallab.com
www.MicroEssentialLab.com

HYDRION PERACETIC ACID (PAA) TEST PAPER, 0-160 AND 0-2000 PPM

- Strip dimensions: 2 1/2 in. x 3/16 in.
- Packaged in flip-top vials
- Color chart label

Catalog No.	pH Range	Color Matching Points
PAA-160	0-160	0-10-30-60-80-160 ppm
PAA-1000	0-1000	0-50-100-250-500-1000 ppm
PAA-2000	0-2000	0-250-500-1000-1500-2000 ppm

Sale Unit	Price	Minimum Order Quantity
Case of 6 Vials	Distributor pricing available	1 case



Scan code to watch instructional videos.

© 2020 Micro Essential Laboratory Inc. All rights reserved. Printed in U.S.A.

For more information contact
Rob Lynch at 718.338.3618
Rob@microessentiallab.com
www.MicroEssentialLab.com

PAA/08.10.20