

Hydrion® Urine & Saliva pH Paper 5.5-8.0

Trusted by healthcare professionals for the fast, easy, and accurate measurement of urine and saliva pH around neutral (pH 7.0).

Widely recommended by physicians and holistic health practitioners for home use (often in conjunction with an alkalizing diet), Hydrion 5.5-8.0 Urine & Saliva pH Paper and Strips set the standard for testing and monitoring body pH. They offer easy-to-read, highly accurate, single color matches at every 0.2 to 0.3 interval from pH 5.5 to 8.0, with excellent color separation in both the acid range (below pH 7.0) and the alkali range (above pH 7.0). Economical and convenient to use for frequent testing.

- Fast, accurate, easy-to-use
- Single color matches
- 9700 Is packaged in a convenient flip-top vials
- 067 comes in a convenient test paper dispenser
- Private labeling options available
- Made in the U.S.A..



Cat# 9700
5.5-8.0 plastic strips



Cat# 067
5.5-8.0 pH paper

For more information contact
Rob Lynch at 718.338.3618
Rob@microessentiallab.com
www.MicroEssentialLab.com



Scan code to watch instructional videos.

HYDRION® (O67) URINE & SALIVA PH PAPER 5.5-8.0

- Each kit consists of a 15 ft. x 1/4 in. roll of premium test paper in a handy Hydrion dispenser, plus a matching color chart
- Approximately 100 tests per kit



catalog no.	pH range	pH readings
O67	5.5-8.0	5.5-5.8-6.0-6.2-6.4-6.6-6.8-7.0-7.2-7.4-7.6-8.0

sale unit	price	minimum order quantity
Case of 10 kits	Distributor pricing available	1 case

HYDRION® (9700) URINE & SALIVA PH PLASTIC STRIPS 5.5-8.0

- Each vial contains 100 strips



catalog no.	pH range	pH readings
9700	5.5-8.0	5.5-5.8-6.0-6.2-6.4-6.6-6.8-7.0-7.2-7.4-7.6-8.0

sale unit	price	minimum order quantity
Case of 6 vials	Distributor pricing available	1 case

PREFERRED BY HEALTHCARE PRACTITIONERS

Many homeopathic nutritionists stress the importance of maintaining an optimum body pH range. And for the regular monitoring of body pH, no test paper is more trusted by physicians and healthcare practitioners than Hydrion (O67) Urine & Saliva pH Paper. Millions of kits sold worldwide attest to the paper's quality, excellent fidelity, precise measurement, and ease of use.

NOTE: Hydrion products are not intended for the diagnosis or treatment of any disease or condition. Contact a licensed physician or healthcare practitioner before drawing conclusions about your health based on the use of urine or saliva pH tests.

© 2019 Micro Essential Laboratory Inc. All rights reserved. Printed in U.S.A.

For more information contact
 Rob Lynch at 718.338.3618
 Rob@microessentiallab.com
 www.MicroEssentialLab.com

5.5-8.0/01.08.19