

## HYDRION® (9800) PLASTIC PH STRIPS 0-14

Used worldwide, Hydrion Plastic pH Strips, set the standard for quick, easy, and accurate testing across the entire pH range.

The latest addition to our bestselling line of Hydrion brand pH test products, Hydrion 0-14 Plastic pH Strips provide excellent color separation over the entire range of the pH scale, with clear, bright single color matches at every 1.0 interval from pH 0.0 to 14.0. The ultimate in easy and accurate pH testing, the new Hydrion Strips are packaged in convenient flip-top vials with matching color charts.

- **Fast, accurate, easy-to-use:  
one dip, one color match**
- **Packaged in convenient flip-top vials**
- **Made in the U.S.A.**



For more information contact  
Rob Lynch at 718.338.3618  
Rob@microessentiallab.com  
www.MicroEssentialLab.com

## Hydrion® (9800) PLASTIC pH STRIPS 0-14

- Strip dimensions: 2 3/4 in. x 1/4 in.
- Packaged in flip-top vials
- 100 strips per vial

Catalog No.	pH Range	pH Readings
9800	0-14	0.0   1.0   2.0   3.0   4.0   5.0   6.0   7.0   8.0   9.0   10.0   11.0   12.0   13.0   14.0

Sale Unit	Price	Minimum Order Quantity
Case of 6 vials	Distributor pricing available	1 case

### How to Measure pH with Hydrion Strips

- Dip the strip into the test solution
- Remove the strip and immediately compare the resulting color with the matching pH color chart



Scan code to watch instructional videos.



© 2019 Micro Essential Laboratory Inc. All rights reserved. Printed in U.S.A.

For more information contact  
 Rob Lynch at 718.338.3618  
 Rob@microessentiallab.com  
 www.MicroEssentialLab.com

ME9800/01.08.19