



**HydriON® Buffer Certification Sheet**  
**&**  
**Certificate of Traceability and Conformity**

**pH Value** 4.00  
**Certified at 25°C** 4.00±0.02

**Mfg. Lot #** 229123 **Exp Date** 10-15-2028

Traceable IUPAC/NIST standard and/or Certified Reference Materials used for Calibration of pH Meters

Buffer pH	NIST Standard Reference Material Used	NIST Standard Reference Material No.
7.00	POTASSIUM DIHYDROGEN PHOSPHATE, DISODIUM HYDROGEN PHOSPHATE	186-I-g 186-II-g
4.005	POTASSIUM HYDROGEN PHTHALATE	185i

**pH Values Recorded**

pH Meter	Calibration Date	Test Date	Sample Weight (grams)	pH value
#1	10-18-2023	10-18-2023	1.0 GRAM	3.99
#2	10-18-2023	10-18-2023	1.0 GRAM	4.00

Comments:

This certificate certifies that the above buffer salt here-in described has been analyzed with appropriate IUPAC/National Institute of Standards and Technology Standard, Certified Reference Materials, or internationally recognized traceable standards.

Certified By: Paresh Shah Signature: *Shah PV* Date: 10-18-2023